

High Purity Porous Metal Products for Biotechnology Applications

When it comes to characteristics like strength, precision, efficiency, longevity, or clean-ability, no other material performs more reliably than Mott porous metal. That's why Mott media is designed into dozens of applications in the biotechnology and pharmaceutical industries.



*Mott In-tank Porous Metal Spargers
for efficiency in mass transfer*

Uniform High Precision Porosity

Strictly controlled porosity in media grades 0.2 μ m to 100 μ m.

Depth Filtration

The compaction and sintering process arranges pores randomly, creating a labyrinth of irregular pathways which effectively capture particles that enter the media. Particle capture efficiencies down to 0.003 μ m are obtained.

Corrosion Resistance

Mott porous metal products are available in corrosion resistant alloys like 316LSS, Nickel, Hastelloy[®], Titanium and other alloys offering long term resistance to corrosive cleaning agents.

Hastelloy is a registered trademark of Haynes International, Inc.

mott corporation

84 Spring Lane, Farmington, CT 06032-3159
www.mottcorp.com, email: quest@mottcorp.com
860-747-6333 Fax 860-747-6739

Tolerance of High Temperatures and Pressures

Mott porous metal products withstand temperatures up to 1700°F (925°C) in reducing or neutral environments and 1450°F (787°C) in oxidizing atmospheres and differential pressures in excess of 3000 psig (206 bar).

Clean-ability

Mott porous metal products can be cleaned in-situ with blowback (gas) or back pulse (liquids) methods. They can be washed or purged with water, solvents, bleach, caustic or acid solutions. They also can be autoclaved, steam purged or cleaned ultrasonically.

Mott porous metal applications in the biotechnology and pharmaceutical industries include:

- Bioreactor/fermentor process gas spargers
- Vent filters
- Steam filters
- Steam spargers
- Gas line filters



*Porous Metal Filter Cartridges for Vent,
Steam and Gas Filtration*



ISO 9001:2000 CERTIFIED

HPBIOAPP REV 1 0109